

Product datasheet for **TA385599**

BrdU Mouse Monoclonal Antibody [Clone ID: 41B2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	41B2
Applications:	IF
Recommended Dilution:	ICC/IF: 1/50-1/200
Reactivity:	All Species
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	Mouse IgG in phosphate buffered saline , pH 7.4, 150mM NaCl and 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store at 2-8 °C. Avoid freeze/thaw cycles.
Stability:	1 year
Background:	BrdU is a uridine derivative and a structural analog of thymidine, which can be incorporated into DNA during the S-phase of a cell cycle as a substitute for thymidine.



[View online »](#)