

Product datasheet for **TA385611**

Recombinant Anti-c-myc epitope tag (Clone 9E10)

Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-c-myc epitope tag (Clone 9E10)
Clone Name:	9E10
Applications:	IF, IHC, IP, WB
Reactivity:	Myc tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
Isotype:	IgG2a, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	v-myc avian myelocytomatosis viral oncogene homolog
Note:	This is a reformatted mouse IgG2a antibody, based on the original mouse IgG1 format, created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)