

## Product datasheet for **TA385620**

### Recombinant Anti-His tag (C-term) (Clone 3D5)


#### Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-His tag (C-term) (Clone 3D5)
Clone Name:	3D5
Applications:	ELISA, WB
Reactivity:	His tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Purified his tagged Fab fragment.
Isotype:	IgG1, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a reformatted mouse IgG1 antibody, based on the original mouse IgG2b format, created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)

## Product images:



TA385620 (1  $\mu\text{g}/\text{ml}$ , right panel) staining of Multitag protein at 25 ng (A), 10 ng (B) and 5 ng (C) amounts. Primary incubation was 1 hour. Detected by chemiluminescence.