

Product datasheet for **TA385630**

Recombinant Anti-His tag (C-term) (Clone 3D5)

Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-His tag (C-term) (Clone 3D5)
Clone Name:	3D5
Applications:	ELISA, WB
Reactivity:	His tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Purified his tagged Fab fragment.
Isotype:	IgG3, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a reformatted mouse IgG3 antibody, based on the original mouse IgG2b format, created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)