

Product datasheet for **TA385635**

FITC Mouse Monoclonal Antibody [Clone ID: 4-4-20 (enhanced)]

Product data:

| | |
|----------------------|---|
| Product Type: | Primary Antibodies |
| Clone Name: | 4-4-20 (enhanced) |
| Host: | Mouse |
| Isotype: | IgG1, kappa |
| Clonality: | Monoclonal |
| Immunogen: | Fluorescein I-bovine serum albumin (FII-BSA). |
| Specificity: | Fluorescein labelled proteins. |

This antibody is not conjugated to fluorescein, but was raised against fluorescein. This antibody can be used as an isotype control in both in vitro and in vivo applications.

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.

Stability: 3 years from dispatch.

Note: This is a reformatted mouse IgG1 antibody, based on the original mouse IgG2a format, created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)