

Product datasheet for **TA385645**

FITC Mouse Monoclonal Antibody [Clone ID: 4-4-20 (enhanced)]

Product data:

Product Type:	Primary Antibodies
Clone Name:	4-4-20 (enhanced)
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Immunogen:	Fluorescein I-bovine serum albumin (FII-BSA).
Specificity:	Fluorescein labelled proteins.

This antibody is not conjugated to fluorescein, but was raised against fluorescein. This antibody can be used as an isotype control in both in vitro and in vivo applications.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a chimeric antibody created as part of a panel offering antibodies of the same specificity in different formats (species, isotype, subtype and modified versions) for use as isotype controls.



[View online »](#)