

## Product datasheet for **TA385650**

### **FITC Rat Monoclonal Antibody [Clone ID: 4-4-20 (enhanced)]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	4-4-20 (enhanced)
<b>Host:</b>	Rat
<b>Isotype:</b>	IgG2b, kappa
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Fluorescein I-bovine serum albumin (FII-BSA).
<b>Specificity:</b>	Fluorescein labelled proteins.

This antibody is not conjugated to fluorescein, but was raised against fluorescein. This antibody can be used as an isotype control in both in vitro and in vivo applications.

**Formulation:** PBS with 0.02% Proclin 300.

**Concentration:** lot specific

**Conjugation:** Unconjugated

**Storage:** Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.

**Stability:** 3 years from dispatch.

**Note:** This is a chimeric antibody created as part of a panel offering antibodies of the same specificity in different formats (species, isotype, subtype and modified versions) for use as isotype controls.



[View online »](#)