

## Product datasheet for **TA385656**

### Hapten 4-hydroxy-3-nitrophenyl acetyl (NP) Mouse Monoclonal Antibody [Clone ID: B1-8]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B1-8
Host:	Mouse
Isotype:	Fab fragment, lambda
Clonality:	Monoclonal
Immunogen:	hapten-protein conjugated NP12-chicken gammaglobulin (CG).
Specificity:	Hapten (4-hydroxy-3-nitrophenyl) acetyl (NP).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	His Tag
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a chimeric antibody created as part of a panel offering antibodies of the same specificity in different formats (species, isotype, subtype and modified versions) for use as isotype controls.



[View online »](#)