

Product datasheet for **TA385673**

Cd3e Monoclonal Antibody [Clone ID: 145-2C11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	145-2C11
Applications:	FC, IP, WB
Reactivity:	Mouse
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	H-2Kb-specific murine cytotoxic T-lymphocyte (CTL) clone.
Specificity:	Recognises a 25-kDa protein component (CD3e, originally called T3-e) of the antigen-specific T-cell receptor.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD3 antigen, epsilon polypeptide
Database Link:	Entrez Gene 12501 Mouse P22646
Synonyms:	CD3-epsilon; FLJ18683; T3E; TCRE



[View online »](#)