

Product datasheet for **TA385679**

CD200R (CD200R1) Mouse Monoclonal Antibody [Clone ID: OX108]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OX108
Applications:	FC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Recombinant protein - hCD200RCD4d3+4.
Specificity:	Recognizes human CD200R, a receptor for the CD200/OX2.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD200 receptor 1
Database Link:	Entrez Gene 131450 Human Q8TD46
Synonyms:	CD200R; CRTR2; HCRTR2; MOX2R; OX2R



[View online »](#)