

Product datasheet for **TA385686**

CD3E Rabbit Monoclonal Antibody [Clone ID: OKT-3 (muromonab)]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OKT-3 (muromonab)
Applications:	FC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	Human peripheral E+ cells.
Specificity:	Epitope on the epsilon-subunit within the human CD3 complex.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD3e molecule
Database Link:	Entrez Gene 916 Human P07766
Synonyms:	CD3-epsilon; FLJ18683; T3E; TCRE
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG2a format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)