

## Product datasheet for **TA385705**

### EGFR Mouse Monoclonal Antibody [Clone ID: 528]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	528
Applications:	FC, IF, IHC, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Partially purified EGFR from A431 cells.
Specificity:	Recognizes the extracellular domain of EGFR.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	epidermal growth factor receptor
Database Link:	<a href="#">Entrez Gene 1956 Human P00533</a>
Synonyms:	ERBB; ERBB1; HER1; mENA; PIG61



[View online »](#)