

Product datasheet for **TA385707**

GLB1 Mouse Monoclonal Antibody [Clone ID: 13R4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	13R4
Applications:	FC, WB
Reactivity:	Escherichia coli
Host:	Mouse
Isotype:	IgG1, lambda
Clonality:	Monoclonal
Immunogen:	Beta-galactosidase.
Specificity:	13R4 recognizes wild-type E. coli beta-galactosidase.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	galactosidase beta 1
Database Link:	P16278
Synonyms:	EBP; ELNR1; Lactase



[View online »](#)