

## Product datasheet for **TA385720**

### GLB1 Rat Monoclonal Antibody [Clone ID: 13R4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	13R4
Applications:	FC, WB
Reactivity:	Escherichia coli
Host:	Rat
Isotype:	IgG1, lambda
Clonality:	Monoclonal
Immunogen:	Beta-galactosidase.
Specificity:	13R4 recognizes wild-type E. coli beta-galactosidase.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	galactosidase beta 1
Database Link:	<a href="#">P16278</a>
Synonyms:	EBP; ELNR1; Lactase
Note:	This is a chimeric antibody created as part of a panel offering antibodies of the same specificity in different formats (species, isotype, subtype and modified versions) for use as isotype controls.



[View online »](#)