

## Product datasheet for **TA385724**

### Recombinant Anti-V5 epitope tag (Clone SV5-P-K)


#### Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-V5 epitope tag (Clone SV5-P-K)
Clone Name:	SV5-P-K
Applications:	ELISA, IF, IHC, IP, R, WB
Reactivity:	V5 tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Pk (V5) Epitope Tag (GKPIPPELLGLDST).
Isotype:	IgG1, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.



[View online »](#)

**Product images:**



ELISA of anti-epitope tag antibodies on multi-tag protein Plate coated with 5ug/ml of Multitag protein. Primary antibody loaded at 1000 ng/ml in first well and an 11 point serial dilution performed in duplicate for each antibody with last well left blank. Signal detection was performed using HRP-conjugated goat anti-mouse.