

Product datasheet for TA385728

Recombinant Anti-V5 epitope tag (Clone SV5-P-K)

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-V5 epitope tag (Clone SV5-P-K)
Clone Name:	SV5-P-K
Applications:	ELISA, IF, IHC, IP, R, WB
Reactivity:	V5 tag
Host:	Human
Clonality:	Monoclonal
Immunogen:	Pk (V5) Epitope Tag (GKPIPNPLLGLDST).
lsotype:	lgG2, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a chimeric human antibody, based on the original mouse IgG1 format, created for improved compatibility with existing reagents, assays and techniques.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US