

Product datasheet for TA385730

Recombinant Anti-V5 epitope tag (Clone SV5-P-K)

Product data:

Product Type: Tag Antibodies

Product Name: Recombinant Anti-V5 epitope tag (Clone SV5-P-K)

Clone Name: SV5-P-K

Applications: ELISA, IF, IHC, IP, R, WB

Reactivity: V5 tag
Host: Human

Clonality: Monoclonal

Immunogen: Pk (V5) Epitope Tag (GKPIPNPLLGLDST).

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid

freeze and thaw cycles.

Stability: 3 years from dispatch.

Note: This is a chimeric human antibody, based on the original mouse IgG1 format, created for

improved compatibility with existing reagents, assays and techniques.

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com