

Product datasheet for **TA385731**

Recombinant Anti-V5 epitope tag (Clone SV5-P-K)

Product data:

Product Type:	Tag Antibodies
Product Name:	Recombinant Anti-V5 epitope tag (Clone SV5-P-K)
Clone Name:	SV5-P-K
Applications:	ELISA, IF, IHC, IP, R, WB
Reactivity:	V5 tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Pk (V5) Epitope Tag (GKPIPPELLGLDST).
Isotype:	IgG2a, kappa
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a reformatted mouse IgG2a antibody, based on the original mouse IgG1 format, created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)