

## Product datasheet for **TA385745**

### IL18 Mouse Monoclonal Antibody [Clone ID: 125-2H]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	125-2H
Applications:	ELISA, IF, IP
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Fusion of splenocytes from mice immunized with human IL-18, with SP2/O mouse myeloma cells.
Specificity:	Recognizes the active form of human IL-18.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	interleukin 18
Database Link:	<a href="#">Entrez Gene 3606 Human Q14116</a>
Synonyms:	lbctadekin; IGIF; IL-1g; IL-18; IL1F4; interleukin-18; MGC12320



[View online »](#)