

Product datasheet for **TA385747**

p53 (TP53) Mouse Monoclonal Antibody [Clone ID: PAb421]

Product data:

| | |
|----------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | PAb421 |
| Applications: | ELISA, IF, IHC, IP, WB |
| Reactivity: | Human, Monkey, Mouse, Rabbit, Rat |
| Host: | Mouse |
| Isotype: | IgG2a, kappa |
| Clonality: | Monoclonal |
| Immunogen: | Synthetic peptide corresponding to aa 371-380 of human p53. |
| Specificity: | Recognises the mammalian mutant and wild-type p53 protein and reacts with an epitope between amino acid residues 370 and 378 (near the C-terminus of the protein), corresponding to the sequence "KKGQSTSRHK". |
| Formulation: | PBS with 0.02% Proclin 300. |
| Concentration: | lot specific |
| Conjugation: | Unconjugated |
| Storage: | Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles. |
| Stability: | 3 years from dispatch. |
| Gene Name: | tumor protein p53 |
| Database Link: | Entrez Gene 7157 Human P04637 |
| Synonyms: | FLJ92943; LFS1; p53; TRP53 |



[View online »](#)