

Product datasheet for **TA385749**

P Glycoprotein (ABCB1) Mouse Monoclonal Antibody [Clone ID: MRK16]

Product data:

Product Type:	Primary Antibodies
Clone Name:	MRK16
Applications:	FC, IF, IP
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Drug resistant K-562 cells.
Specificity:	Recognizes a discontinuous extracellular epitope on P-glycoprotein.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	ATP binding cassette subfamily B member 1
Database Link:	Entrez Gene 5243 Human P08183
Synonyms:	ABC20; CD243; CLCS; GP170; IBD13; MDR1; MGC163296; P-GP; PGY1



[View online »](#)