

## Product datasheet for **TA385761**

### CD2 Rat Monoclonal Antibody [Clone ID: YTH 655]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	YTH 655
Applications:	FC
Reactivity:	Human
Host:	Rat
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Specificity:	This antibody specifically recognizes human CD2, a cell adhesion molecule found on the surface of T cells and natural killer (NK) cells.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD2 molecule
Database Link:	<a href="#">Entrez Gene 914 Human P06729</a>
Synonyms:	FLJ46032; LFA-2; SRBC; T11



[View online »](#)