

Product datasheet for **TA385774**

CD3E Rat Monoclonal Antibody [Clone ID: YTH 12.5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	YTH 12.5
Applications:	IF, IHC, WB
Reactivity:	Human
Host:	Rat
Isotype:	IgG2b, lambda
Clonality:	Monoclonal
Immunogen:	Human CD3.
Specificity:	Recognises the epsilon chain of the human CD3 antigen complex on human T cells.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD3e molecule
Database Link:	Entrez Gene 916 Human P07766
Synonyms:	CD3-epsilon; FLJ18683; T3E; TCRE



[View online »](#)