

## Product datasheet for **TA385779**

### NOTCH2 Mouse Monoclonal Antibody [Clone ID: B6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B6
Applications:	BI, ELISA, IF
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	NRR of mouse Notch 2.
Specificity:	An engineered ScFv fragment antibody that binds to both mouse mNRR2 and human hNRR2, but not NRR1.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	notch 2
Database Link:	<a href="#">Entrez Gene 4853 Human Q04721</a>
Synonyms:	AGS2; hN2
Note:	This is a chimeric antibody created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)