

Product datasheet for **TA385783**

Recombinant Anti-Unknown Specificity (Isotype control) (Clone MOPC-21) Isotype Control

Product data:

Product Type:	Isotype Controls
Clone Name:	MOPC-21
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Specificity:	An antibody of unknown specificity that maybe used as a negative control.
Buffer:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.



[View online »](#)