

## Product datasheet for **TA385792**

### Recombinant Anti-Unknown Specificity (Isotype control) (Clone MOPC-21) Isotype Control

#### Product data:

Product Type:	Isotype Controls
Clone Name:	MOPC-21
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Specificity:	An antibody of unknown specificity that maybe used as a negative control.
Buffer:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This is a chimeric antibody created as part of a panel offering antibodies of the same specificity in different formats (species, isotype, subtype and modified versions) for use as isotype controls.



[View online »](#)