

Product datasheet for **TA385804**

beta 2 Microglobulin (B2M) Mouse Monoclonal Antibody [Clone ID: BBM.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	BBM.1
Applications:	IF, IHC, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Immunogen:	Molt 4, a human T cell line.
Specificity:	Binds human beta-2-microglobulin, a component of MHC class I molecules, which are present on all nucleated cells (excludes red blood cells) and involved in the presentation of peptide antigens to the immune system.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	beta-2-microglobulin
Database Link:	Entrez Gene 567 Human P61769
Synonyms:	beta-2-microglobin; beta-2-microglobulin



[View online »](#)