

Product datasheet for **TA385807**

HLA class II Mouse Monoclonal Antibody [Clone ID: F3.3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	F3.3
Applications:	FC, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Specificity:	Recognising anti-class II HLA molecules.

This antibody has been shown to work by flow-cytometry and western blot. Dechant et al (PMID: 12393717) used a human IgA2-chimeric version of this antibody to study IgA-mediated cell killing.

Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.



[View online »](#)