

Product datasheet for **TA385815**

alpha 1 Fetoprotein (AFP) Mouse Monoclonal Antibody [Clone ID: AF5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AF5
Applications:	FC, IF, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	AF5 was prepared by immunizing female BALB/c mice with human amniotic fluid containing AFP, from a pregnancy affected by a neural tube defect, diluted and emulsified in Freund's complete adjuvant.
Specificity:	Recognises the human Alpha-fetoprotein.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	alpha fetoprotein
Database Link:	Entrez Gene 174 Human P02771
Synonyms:	alpha-1-fetoprotein; Alpha-fetoglobulin; alpha-fetoprotein; FETA; HP; HPAFP; OTTHUMP00000160480



[View online »](#)