

## Product datasheet for **TA385817**

### CD19 Mouse Monoclonal Antibody [Clone ID: HD37]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	HD37
Applications:	Bl, FC, IF, IHC, IP
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Cells from a hairy cell leukemia patient.
Specificity:	Recognises human CD19.
	This antibody binds to CD19 and is able to inhibit B-cell activation (Pezutto 1987).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD19 molecule
Database Link:	<a href="#">Entrez Gene 930 Human P15391</a>
Synonyms:	B4; Leu-12; MGC12802



[View online »](#)