

Product datasheet for **TA385823**

CD19 Mouse Monoclonal Antibody [Clone ID: 4G7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	4G7
Applications:	FC, IF
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Human CCL (chronic lymphocytic leukemia) cells.
Specificity:	Recognises human CD19.
	This antibody was used in the identification of the B-cell epitope CD19.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD19 molecule
Database Link:	Entrez Gene 930 Human P15391
Synonyms:	B4; Leu-12; MGC12802



[View online »](#)