

## Product datasheet for **TA385825**

### Integrin beta 1 (ITGB1) Mouse Monoclonal Antibody [Clone ID: K20]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	K20
Applications:	FC, IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Human T-cell lymphoma cells injected into Biozzi high-responder mice.
Specificity:	Recognises human CD29.  This antibody binds to CD29.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	integrin subunit beta 1
Database Link:	<a href="#">Entrez Gene 3688 Human P05556</a>
Synonyms:	CD29; FNRB; GPIIA; MDF2; MSK12; VLA-BETA; VLAB



[View online »](#)