

## Product datasheet for **TA385828**

### CCR3 Mouse Monoclonal Antibody [Clone ID: 5E8-G9-B4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	5E8-G9-B4
Applications:	FC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Specificity:	Recognises human CCR3, a chemokine receptor highly expressed of eosinophils, as well as Th2-Tcells, basophils, platelets, mast-cells, mononuclear phagocytes, airway epithelial cells, and hematopoietic progenitor cells (CD34+).
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	C-C motif chemokine receptor 3
Database Link:	<a href="#">Entrez Gene 1232 Human P51677</a>
Synonyms:	CC-CKR-3; CCR-3; CD193; CKR3; CMKBR3; MGC102841; OTTHUMP00000164718; OTTHUMP00000209748



[View online »](#)