

## Product datasheet for **TA385831**

### IL2 Mouse Monoclonal Antibody [Clone ID: LNKB-2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	LNKB-2
Applications:	ELISA, IF, IP, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Specificity:	Recognises human Interleukin 2.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	interleukin 2
Database Link:	<a href="#">Entrez Gene 3558 Human P60568</a>
Synonyms:	aldesleukin; IL-2; interleukin-2; lymphokine; TCGF



[View online »](#)