

Product datasheet for **TA385835**

CD37 Mouse Monoclonal Antibody [Clone ID: WR17]

Product data:

Product Type:	Primary Antibodies
Clone Name:	WR17
Applications:	FC, IF, IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Chronic lymphocytic leukaemia cells.
Specificity:	Recognises human CD37.
	This antibody binds CD37, staining, amongst others, B-cells and B-cell tumours.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD37 molecule
Database Link:	Entrez Gene 951 Human P11049
Synonyms:	GP52-40; MGC120234; tetraspanin-26; Tspan-26; TSPAN26



[View online »](#)