

Product datasheet for **TA385842**

CD4 Mouse Monoclonal Antibody [Clone ID: 13B8.2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	13B8.2
Applications:	EIA, ELISA, FC, IHC, WB
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Human HSB2-T cells.
Specificity:	Recognises human CD4.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD4 molecule
Database Link:	Entrez Gene 920 Human P01730
Synonyms:	CD4mut



[View online »](#)