

## Product datasheet for **TA385845**

### CD7 Mouse Monoclonal Antibody [Clone ID: 3A1E]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	3A1E
Applications:	FC, Func
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Specificity:	Recognises human CD7.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD7 molecule
Database Link:	<a href="#">Entrez Gene 924 Human P09564</a>
Synonyms:	GP40; Leu-9; Tp40; TP41



[View online »](#)