

Product datasheet for **TA385848**

EGFR Rabbit Monoclonal Antibody [Clone ID: C225 (Cetuximab)]

Product data:

Product Type:	Primary Antibodies
Clone Name:	C225 (Cetuximab)
Applications:	BI, FC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	The epidermal growth factor receptor is the cell-surface receptor for members of the epidermal growth factor family (EGF-family) of extracellular protein ligands.
Specificity:	Recognises human EGFR.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	epidermal growth factor receptor
Database Link:	Entrez Gene 1956 Human P00533
Synonyms:	ERBB; ERBB1; HER1; mENA; PIG61
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original human IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)