

## Product datasheet for **TA385852**

### CEA Rabbit Monoclonal Antibody [Clone ID: Arcitumomab]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	Arcitumomab
Applications:	FC, IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	Carcinoembryonic antigen (CEA) describes a set of highly related glycoproteins involved in cell adhesion normally produced in gastrointestinal tissue during fetal development, but only at very low levels in adults.
Specificity:	Recognises human carcinoembryonic antigen.  This antibody binds to carcinoembryonic antigen and was used in the diagnosis of mammary carcinomas.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Database Link:	<a href="#">Q13984</a>
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original human IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)