

Product datasheet for **TA385863**

Fcgr2b Mouse Monoclonal Antibody [Clone ID: AT130-5]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AT130-5
Applications:	IHC
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	CD32-rat CD4 domain 3/4 fusion protein.
Specificity:	Binds specifically to mouse CD32. The antibody binds to CD32, a surface receptor protein with a low affinity for IgG antibodies.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	Fc receptor, IgG, low affinity IIb
Database Link:	Entrez Gene 14130 Mouse P08101
Synonyms:	CD32; CD32B; CDw32; Fc-gamma-RIIb; FCG2; FCGR2; FcRII-b; IGFR2



[View online »](#)