

Product datasheet for **TA385865**

CD32B (FCGR2B) Mouse Monoclonal Antibody [Clone ID: AT10]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AT10
Applications:	FC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	Human CD32.
Specificity:	Binds specifically to human CD32.
	The antibody binds to CD32, a surface receptor protein with a low affinity for IgG antibodies.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	Fc fragment of IgG receptor IIb
Database Link:	Entrez Gene 2213 Human P31994
Synonyms:	CD32; CD32B; CDw32; Fc-gamma-RIIb; FCG2; FCGR2; FcRII-b; IGFR2



[View online »](#)