

## Product datasheet for **TA385871**

### IgM Mouse Monoclonal Antibody [Clone ID: M15/8]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	M15/8
Applications:	BI, ELISA, FC, IHC, IP
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	BALB/c mice immunized with human IgM.
Specificity:	This antibody binds to human IgM.

This antibody binds to IgM, which acts in primary immune defines and is involved in early recognition of antigens.

**Formulation:** PBS with 0.02% Proclin 300.

**Concentration:** lot specific

**Conjugation:** Unconjugated

**Storage:** Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.

**Stability:** 3 years from dispatch.

**Database Link:** [P01871](#)



[View online »](#)