

## Product datasheet for **TA385876**

### Syntaxin Rabbit Monoclonal Antibody [Clone ID: SP6]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	SP6
Applications:	IHC, WB
Reactivity:	Human, Rat, Hamster, Rabbit, Pig, Dog, Guinea Pig, Cow
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	Synaptic vesicle-containing fractions immunoprecipitated from human brain homogenates.
Specificity:	This antibody binds to native and denatured syntaxin.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)