

Product datasheet for **TA385878**

SNAP25 Rabbit Monoclonal Antibody [Clone ID: SP12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	SP12
Applications:	ELISA, IF, IHC, WB
Reactivity:	Human, Mouse, Rat, Hamster, Pig, Feline, Guinea Pig, Cow
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Specificity:	This antibody binds to native and denatured SNAP-25.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	synaptosome associated protein 25kDa
Database Link:	Entrez Gene 6616 Human P60880
Synonyms:	bA416N4.2; dj1068F16.2; FLJ23079; RIC-4; RIC4; SEC9; SNAP; SNAP-25; SUP
Note:	This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG1 format, for improved compatibility with existing reagents, assays and techniques.



[View online »](#)