

Product datasheet for TA385888

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GRB2 Mouse Monoclonal Antibody [Clone ID: SAIC-19D-1F4]

Product data:

Product Type: Primary Antibodies

Clone Name: SAIC-19D-1F4

Applications: WB

Reactivity: Human
Host: Mouse

Isotype: IgG1, kappa
Clonality: Monoclonal

Immunogen: Peptide "FGNDVQHFK", derived from the Growth factor receptor-bound protein 2, conjugated

to KLH.

Specificity: Recognises human GRB2.

Original data characterizing this antibody may be foundhere.

Formulation: PBS with 0.02% Proclin 300.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid

freeze and thaw cycles.

Stability: 3 years from dispatch.

Gene Name: growth factor receptor bound protein 2

Database Link: Entrez Gene 2885 Human

P62993

Synonyms: AA408164; Ash; MGC108668

Note: This is a chimeric antibody created for improved compatibility with existing reagents, assays

and techniques.

