

Product datasheet for **TA385894**

IL18 Mouse Monoclonal Antibody [Clone ID: SAIC-24B-4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	SAIC-24B-4
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1, kappa
Clonality:	Monoclonal
Immunogen:	The KLH conjugated peptide "ISTLSCENK" derived from IL-18, an augmentor of natural killer cell activity in the spleen and stimulator of interferon gamma production in Th1 cells.
Specificity:	Recognises human IL18. Original data characterizing this antibody may be found here .
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	interleukin 18
Database Link:	Entrez Gene 3606 Human Q14116
Synonyms:	lboctadekin; IGIF; IL-1g; IL-18; IL1F4; interleukin-18; MGC12320
Note:	This is a chimeric antibody created for improved compatibility with existing reagents, assays and techniques.



[View online »](#)