

## Product datasheet for **TA385898**

### **NFkB p100 / p52 (NFKB2) Rabbit Monoclonal Antibody [Clone ID: SAIC-26C-15]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	SAIC-26C-15
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	Peptide "IEVDLVTHSDPPR" derived from the transcription factor NFKB2 (Nuclear factor NF-kappa-B p100 subunit) conjugated to KLH.
Specificity:	Recognises human NFKB2.  Original data characterizing this antibody may be found <a href="#">here</a> .
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	nuclear factor kappa B subunit 2
Database Link:	<a href="#">Entrez Gene 4791 Human Q00653</a>
Synonyms:	H2TF1; LYT-10; LYT10; OTTHUMP00000020371; p52



[View online »](#)