

Product datasheet for **TA385925**

Nt5e Rat Monoclonal Antibody [Clone ID: YDBYT 208.3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	YDBYT 208.3
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b, kappa
Clonality:	Monoclonal
Immunogen:	Murine CD73.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	5' nucleotidase, ecto
Database Link:	Entrez Gene 23959 Mouse Q61503
Synonyms:	CD73; E5NT; eN; eNT; NT; NT5; NTE



[View online »](#)