

Product datasheet for **TA385927**

PD-L1 (Cd274) Rat Monoclonal Antibody [Clone ID: YDC 127.1.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	YDC 127.1.1
Applications:	FC
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a, kappa
Clonality:	Monoclonal
Immunogen:	Murine PDL1.
Specificity:	This antibody binds to murine PDL1. This antibody binds to murine PDL1.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	CD274 antigen
Database Link:	Entrez Gene 60533 Mouse Q9EP73
Synonyms:	B7-H; B7-H1; B7H1; MGC142294; MGC142296; PD-L1; PDCD1L1; PDCD1LG1; PDL1



[View online »](#)