

## Product datasheet for **TA385929**

### RTN4RL1 Rabbit Monoclonal Antibody [Clone ID: M5]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	M5
Applications:	ELISA, FC, IP, WB
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG, kappa
Clonality:	Monoclonal
Immunogen:	(His)6-tagged fragment of rat NgR1 (amino acids 278-473).
Specificity:	The antibody binds to human and rat NgR1 with a Kd of 2.9 nM for both species (bivalent binding).  The antibody binds specifically to residues 379-448 of NgR1, a member of the Nogo-66 receptor family which inhibits axonal growth and sporulating in the CNS.
Formulation:	PBS with 0.02% Proclin 300.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Please store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. Avoid freeze and thaw cycles.
Stability:	3 years from dispatch.
Gene Name:	reticulon 4 receptor like 1
Database Link:	<a href="#">Entrez Gene 146760 Human Q86UN2</a>
Synonyms:	DKFZp547J144; NgR3; NGRH2; NGRL2



[View online »](#)